









A1320-1



A1320-2



A1320-3



SPECIFICATIONS:

SHUTTERS FIT BRANDSTROM INSTRUMENTS 7/16" SOLENOIDS
WITH $\phi .1247$ SHAFT WITH FLAT, MINIMUM 45° ROTATION

COUNTERWEIGHT MATERIAL: 303 STAINLESS STEEL
FINISH: PASSIVATE

SHUTTER MATERIAL: 5052-H32 ALUMINUM .016" THICK.
FINISH: MATTE BLACK ANODIZE PER MILL-A-8625F, TYPE II, CLASS II

DIMENSION TOLERANCES UNLESS OTHERWISE SPECIFIED				Brandstrom Instruments		
ANGLES	$\pm 1^{\circ}$	2 PLACES	$\pm .010$			
FRACTIONS	$\pm 1/32$	3 PLACES	$\pm .005$			
1 PLACE	$\pm .015$	4 PLACES	$\pm .0005$			
MATERIAL:		SEE SPECS		DRAWN:		P.T.
FINISH:		SEE SPECS		CHECKED:		P.T.
DATE:	6/14/17	SCALE:	NONE	APPROVED:		T.G.A.
TITLE:		SHUTTER ASSEMBLY		DWG NO:	A1320	SHEET: 1 of 1
						REV:

PROPRIETARY DOCUMENT