




D1950-1
D1950-2




D1949-1
D1949-2








D1948-1
D1948-2



D1947-1
D1947-2



D1946-1
D1946-2



NOTE FOR FINISH:

- 1) NUMBER SUFFIX -1: BLACK ANODIZE PER MIL-A-8625F, CLASS II, TYPE II
- 2) NUMBER SUFFIX -2 : NONE

FOR USE ON:

THESE SHUTTERS ARE DESIGNED FOR BOTH BRANDSTROM INSTRUMENTS 3/8" AND 7/16" SOLENOIDS

(1) 2-56 SET SCREW PER SHUTTER PROVIDED FOR INSTALLATION
IF USING ON BRANDSTROM INSTRUMENTS 3/8" SOLENOIDS YOU WILL NEED TO ADD SHUTTER INSERT PART NUMBER D2002 TO ORDER

PROPRIETARY DOCUMENT

DIMENSION TOLERANCES UNLESS OTHERWISE SPECIFIED				Brandstrom Instruments			
ANGLES ± 1°		2 PLACES ±.010					
FRACTIONS ± 1/32		3 PLACES ±.005					
1 PLACE ± .015		4 PLACES ±.0005					
MATERIAL:		6061-T6 ALUMINUM		DRAWN:		P.T.	
FINISH:		SEE NOTES		CHECKED:		P.T.	
DATE: 6/1/16		SCALE: NONE		APPROVED:		T.G.A.	
TITLE:		ACCESSORY SHUTTERS		DWG NO:		SHEET:	REV:
				BI-A-S		1 of 1	